

Innovative solutions from


BucketBuddy™

Magnetic Bucket and Tube Rack for Laboratory Centrifuges

Part# Z326-BKR



 **Magnetic**



- The **BucketBuddy™** is designed to hold all Hermle buckets/inserts up to 750 ml.
- Rack positions hold most capped/rimmed tube sizes from 1.5 ml to 50 ml.
- Modular and Magnetic, use the **BucketBuddy™** in the place and orientation that works for you!
- Save valuable benchspace and time!

Holds buckets/inserts and tubes

Create the layout that works for you!

