

Z287-A Microcentrifuge PLUS

- Bundle includes 2 rotors
 - 18 x 1.5 / 2.0 ml and 2 x microplate
- Up to **14,000 rpm** and **16,058 xg**
- Indicator Light, easily view cycle status

Bundle Includes



Z287-02MP
2 x microplate
3,800 rpm / 1,453 xg



Z287-1820
18 x 1.5 / 2.0 ml
14,000 rpm / 15,994 xg

Also Available



6 x 50 ml
6,000 rpm / 4,427 xg



24 x 1.5 / 2.0 ml
12,000 rpm / 13,683 xg



Unbeatable value!
\$3795

Compact, only 14in. wide

BUNDLE ORDERING INFORMATION

	Reg.	Bundle
Z287-BND-MB	Z287-A Micro bundle (incl. Z287-1820 and Z287-02MP)	\$5530 \$3795

Benchtop Universal Centrifuge Bundles



As low as **\$6,750**

Z326 Series

- Max. swing-out capacity:
8 x 50ml or **12 x 15ml**
- Truly universal, **17** rotor options

BUNDLE ORDERING INFORMATION

		Reg.	Bundle
Z326-BND	Z326 Ambient Swing Out Bundle	\$7,618	\$6,750
Z326-K-BND	Z326-K Refrigerated Swing Out Bundle	\$10,136	\$9,100



As low as **\$8,200**

Z366 Series

- Max. swing-out capacity:
16 x 50ml or **40 x 15ml**
- High speed rotors up to **24,000 xg**

Z366-BND	Z366 Ambient Swing Out Bundle	\$9,087	\$8,200
Z366-K-BND	Z366-K Refrigerated Swing Out Bundle	\$12,136	\$10,600



As low as **\$9,800**

Z446 Series

- Max. swing-out capacity:
28 x 50ml or **68 x 15ml**
- Microplate capacity:
20 x microplates or **8 x PCR Plates**

Z446-BND	Z446 Ambient Swing Out Bundle	\$10,861	\$9,900
Z446-K-BND	Z446-K Refrigerated Swing Out Bundle	\$14,121	\$13,000
Z446-MPBND	Z446 Ambient Microplate Bundle	\$10,903	\$9,800
Z446-K-MPBND	Z446-K Refrigerated Microplate Bundle	\$14,183	\$12,500



*Additional bundles available upon request.
To order 230V voltage, add "E" to part number.

Z306 Universal Centrifuge



- Includes swing-out rotor for up to **8 x 50ml** or **12 x 15ml**
- Universal, 9 rotor options
- Compact, less than 14 inches wide



BUNDLE ORDERING INFORMATION

Z306-BND	Z306 Universal centrifuge with swing-out rotor	Reg. \$5,877	Bundle \$5,300
----------	--	-------------------------	-----------------------

Z216-MK High Speed Refrigerated Microcentrifuge



- **2** bundle options!
- Compact, only 11 inches wide
- High speed, up to **21,631xg**
- **44** place rotor option
- **COMBI-rotor™** option with Quick-Seal™ lid

NEW Hermle COMBI-rotor™ with Quick-Seal lid



Holds 24 x 1.5/2.0ml and 4-PCR strips

Quick-Seal™ lid
Lid opens and locks with only a quarter turn



Now included in all Z216 series bundles!

BUNDLE ORDERING INFORMATION

Z216-K-BND44	Z216-MK Refrigerated Microcentrifuge w/ 44 place rotor	Reg. \$6,569	Bundle \$5,600
Z216-K-CMB	Z216-MK Refrigerated Microcentrifuge w/ COMBI-Rotor™	Reg. \$7,120	Bundle \$5,800

*Ambient versions also available contact your distributor for pricing

Z36 HK Super Speed Centrifuge

*Additional bundles available upon request. To order 230V voltage, add "-E" to part number.



- Speeds up to **65,395 xg!**
- Max. swing-out capacity: **12 x 50ml** or **28 x 15ml**
- Refrigerated, **-20°C to 40°C**



6 x 250ml rotor
10,000 rpm / 15,650 xg



6 x 50ml rotor
21,000 rpm / 41,410 xg



12 x 1.5/2.0 ml rotor
30,000 rpm / 65,395 xg

Purchase a Z36-HK BUNDLE and receive a **FREE BenchPACK!™**

BUNDLE ORDERING INFORMATION

Z36-BND250	Z36-HK Centrifuge with 6 x 250ml rotor	Reg. \$18,652	Bundle \$15,000
Z36-BND50	Z36-HK Centrifuge with 6 x 50ml rotor	Reg. \$16,923	Bundle \$14,000
Z36-BND20	Z36-HK Centrifuge with 12 x 1.5/2.0ml rotor	Reg. \$16,131	Bundle \$14,000

INCLUDES:

- 1 **myFuge™ 12** Mini Centrifuge
- 1 **BenchMixer™ V2** Vortex Mixer
- 1 **myTemp™ Mini** Incubator
- 1 **Accuris™ Mini** Balance
- 1 **Backpack**

Value \$1,500

