

Offers valid until June. 30th 2025

Z206-A



## Z206-A Compact Bundle

- 6 x 50ml rotor



6,000 rpm / 4,427 xg

Z206-BND: **\$2,600**

Z216-MK



## Z216-M Microcentrifuge Bundles

- 44 x 01.5/2.0 ml
- COMBI-rotor, 24 x 1.5/2.0ml and 4 x PCR
- Over 21,000 xg! (15,000 rpm)



44 x 1.5/2.0ml



24 x 1.5/2.0ml  
and 4-PCR strips

**Ambient**  
Z216-BND44:  
(44-place) **\$3,500**

Z216-CMB:  
(COMBI-rotor) **\$4,000**

**Refrigerated**  
Z216-KBND44:  
(44-place) **\$6,200**

Z216-KCMB:  
(COMBI-rotor) **\$6,700**



Z287-A



## Z287-A MicroBundle

- 2 x microplate and 18 x 1.5/2.0ml




3,800 rpm / 3,717 xg



14,000 rpm / 15,994 xg

Z287-BND-MB: **\$4,800**

Z307



## Z307 Deepwell Plate Bundle

- 2 x Deepwell Plate rotor



4,500 rpm / 2,716 xg

Z307-BND-MP: **\$5,500**

Z327-K



## Z327 Series Bundles

- Tissue culture or blood tube options
- up to 8 x 50ml and 16 x 15ml
- up to 28 x 8-10ml and 32 x 1.8-7ml



Tissue Culture Pkg.



Blood Tube Pkg.

5,000 rpm / 4,136 xg

**Ambient**  
Z327-BND:  
(Tissue Culture) **\$8,200**

Z327-BND-BD:  
(Blood tube) **\$8,200**

**Refrigerated**  
Z327-K-BND:  
(Tissue Culture) **\$11,000**

Z327-K-BND-BD:  
(Blood tube) **\$11,000**

Offers valid until June. 30th 2025

Z366-K



## Z366 Series Bundles

- Ambient or Refrigerated
- up to 16 x 50ml and 40 x 15ml



4,500 rpm / 3,735 xg

**Ambient**  
Z366-BND: **\$9,000**

**Refrigerated**  
Z366-K-BND: **\$11,500**

Z446-K



## Z446 Series Bundles

- Tissue culture or Microplate options
- up to 28 x 50ml and 66 x 15ml
- up to 20 std. plates and 4 x Deepwell



Ambient:  
4,000 rpm / 3,434 xg

Refrigerated:  
4,500 rpm / 4,346 xg

Tissue Culture Pkg.    Microplate Pkg.

**Ambient**  
Z446-BND: **\$10,700**  
(Tissue Culture)  
Z446-MPBND: **\$10,700**  
(Microplate)

**Refrigerated**  
Z446-K-BND: **\$14,000**  
(Tissue Culture)  
Z446-K-MPBND: **\$14,000**  
(Microplate)




Z496-K



## Z496-K Floor Standing Bundle

- Swing out rotor: up to 40 x 50ml and 96 x 15ml
- Refrigerated



4,200 rpm / 3,964 xg

Only available in 230V

Z496-K-BND: **\$26,000**

Z36-HK



## Z36-HK Super Speed Bundles

- Max speed up to 65,395 xg
- 3 rotor options available
- Refrigerated



6 x 250ml rotor  
10,000 rpm / 15,650 xg



6 x 50ml rotor  
21,000 rpm / 41,410 xg



12 x 1.5/2.0 ml rotor  
30,000 rpm / 65,395 xg

Z36-BND250: **\$16,000**  
(6 x 250 ml)

Z36-BND50: **\$15,000**  
(6 x 50 ml)

Z36-BND20: **\$14,000**  
(12 x 1.5/2.0 ml)

Z600-S3



## SIEVA III Filtration Bundle

- Includes perforated basket \*filters sold separately
- One-of-a-kind benchtop filtration centrifuge
- Option with temp. equalization port



Perforated rotor  
10,000 rpm / 7,825 xg

\*Available filter pore sizes: 7µ, 20µ, 30µ, 50µ

Z600-S3-BND: **\$20,000**

Z600-S3-TBND: **\$21,000**  
(w/ temp. equalization)